Table 1: Implemented voice commands and their functionalities

Type	Voice Command	Functionality
2*Patient information	patient history	Activates/Deactivates patient history panel
	CT / CT Image / Tomography	Activates/Deactivates CT Visualizer
8*3D visualization	arteries/veins	Activates/Deactivates rendering of all arteries/veins
	Sternum	8*Activates/Deactivates rendering of the associated structure
	Celiac trunk	
	GDA	
	Mesenteric artery	
	Splenic artery	
	Gastric artery	
	Hepatic artery / liver artery	
	Portal vein	
	Vena cava	
	Splenic vein	
	Mesenteric vein	
	Tumor / lesion / cancer	
	Variation	Activates/Deactivates rendering of unusual anatomical structures
3*Control commands	Go up/down	Activates the automatic scrolling up/down of the CT slices
	Stop	Stops the automatic scrolling of CT slices
	Capture photo / hologram	captures a photo using HoloLens mounted camera with / without holograms